SPECTRUM MEDICAL

Surgical Date:	Rehab Start Date:	Distal Bicep Repair Protocol

Distal Bicep Repair Protocol

PHASE 1: 1-3 WEEKS

CLINICAL GOALS:

- Elbow ROM from 30° of extension to 120° of flexion
- · Maintain minimal swelling and soft tissue healing
- Achieve full forearm supination/pronation

TESTING:

Bilateral elbow and forearm ROM

EXERCISES:

- Six times per day the patient should set his IROM hinged brace at 30° of extension and 120° of flexion and perform active assistive flexion and active extension exercises within the brace. Two sets of 10 are performed.
 - The brace is then reset at 90°, the forearm straps are loosened, and the patient performs 2 sets of 10 forearm rotations. The straps should then be secured.
- Ice after exercise, 3-5 times per day
- A sling should be worn only as needed for comfort with the patient maintaining full shoulder ROM.

PHASE II: 3-6 WEEKS

CLINICAL GOALS:

Full elbow and forearm ROM by 6 weeks ◆ Scar management

TESTING:

- Bilateral elbow and forearm ROM
- Grip strengthening at 5-6 weeks

EXERCISES

3 weeks:

- The extension limit in the brace is changed to 20°. Flexion remains at 120°, but patient may remove brace to allow full flexion 2 times per day. The brace stays on at all other times except when washing the arm.
- Scar massage 3-4 times per day.

4 weeks:

- Extension limit is changed to 10°.
- Continue the same exercises.
- Putty may be used 3 times per day for 10 minutes to improve grip strength.

5 weeks:

• The extension limit is changed to 0° and exercises are continued in the brace.

6 weeks:

- The brace is discontinued, unless needed for protection.
- Passive elbow extension exercises are initiated if needed.
- Light strengthening exercises are initiated with light tubing or 2-3 pound weights for elbow flexion, extension, forearm rotation and wrist flexion and extension.
- Ice is continued after strengthening exercises.

CLINICAL FOLLOW-UP:

• The patient usually is seen at 3 weeks and at 5-6 weeks, then only as needed with doctor appointments, until full motion is achieved and to monitor the patient's strengthening program.

SPECTRUM MEDICAL

Surgical Date:	Rehab Start Date:	Distal Bicep Repair Protocol	
PHASE III: 6 WEEKS TO 6 MONTHS			

CLINICAL GOALS:

The strengthening program is gradually increased so that the patient is using full
weights by 3 months. It may be as long as 6 months before a patient returns to
heavy work.

TESTING:

- Grip strengthening
- Elbow ROM

EXERCISES:

- Elbow ROM exercises are performed if ROM is not WNL
- Strengthening exercises to wrist, elbow, forearm, and possibly shoulder depending on sport and/or work requirements

CLINICAL FOLLOW-UP:

 The patient is seen only as needed, usually with doctor appointments, to monitor progress with strengthening program.